

Unistat® 705w

Cooling a Buchi Glas Uster 1-litre reactor to T_{min}

Requirement

This case study looks at the performance of a Unistat 705w cooling a Buchi Glas Uster 1-litre un-insulated jacketed glass pressure reactor to T_{min} under "process" control.

Method

The Unistat and reactor are connected using two 1-metre insulated metal hoses. The reactor is filled with 0.75 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The jacket cooling curve is almost linear to $-50\text{ }^{\circ}\text{C}$ where it begins to asymptote before reaching its lowest temperature of $-66\text{ }^{\circ}\text{C}$ with a corresponding process temperature of $-65\text{ }^{\circ}\text{C}$.

Setup details

Unistat® 705w & Buchi Glas Uster reactor

Temperature range: $-75\text{...}250\text{ }^{\circ}\text{C}$
 Cooling power: $0.6\text{ kW @ }250\text{...}100\text{ }^{\circ}\text{C}$
 $0.65\text{ kW @ }0\text{ }^{\circ}\text{C}$
 $0.6\text{ kW @ }-20\text{...}-40\text{ }^{\circ}\text{C}$
 $0.3\text{ kW @ }-60\text{ }^{\circ}\text{C}$

Heating power: $1.5\text{ kW}/3\text{ kW}$
 Pump speed: 3300 rpm
 Hoses: $2 \times 1\text{ m; M24} \times 1.5\text{ (#9325)}$
 HTF: DW-Therm (#6479)
 Reactor: 1-litre un-insulated jacketed glass pressure reactor

Reactor content: $0.75\text{ litre M90.055.03 (#6259)}$

Stirrer speed: 500 rpm
 Control: process

