

Unistat® 705w

Controlling a Heidolph "Synthesis 1"

Requirement


This test demonstrates the responsiveness and level of control that can be achieved when a Unistat 705 is used in conjunction with a Heidolph "Synthesis 1" parallel synthesis apparatus.

Method

The setpoints of all four temperature zones of the Synthesis 1 were set at -30 °C and then reduced to -50 °C. The response was recorded using the Huber "SpyControl".

Results

The first set-point of -30 °C was reached in approximately 32 min. The second set-point of -50 °C was reached after an additional 23 min. At both set-points the process temperature offset was due to the positioning of the process-located Pt 100.



Unistat® 705w

In comparison to a previously published case study with the Heidolph Synthesis 1 it is possible to achieve lower temperature with less cooling power and in a much faster time demonstrating more efficient thermal transfer.

Setup details

Unistat® 705w with optional bypass installed

- Temperature range: -75...250 °C
- Cooling power: 0.6 kW @ 250...100 °C
0.65 kW @ 0 °C
0.6 kW @ -20...-40 °C
0.3 kW @ -60 °C
- Heating power: 1.5 kW (3 kW @ 400 V)
- Pump speed: 3500 rpm
- Hoses: 2 x 1 m M24 x 1,5 #9325
- HTF: Ethanol
- Reactor: Heidolph Synthesis 1
- Control: process (inside test tube)
- Control of the Heidolph Synthesis 1: off



Setup details: Heidolph Synthesis 1

Each of the 24 test tubes of the Synthesis 1 were filled with 5 ml Sil Oil Huber M90.055.03 #6259

- Control of the Heidolph Synthesis 1: Zone temperatures controlled via internal sensors.
- Setpoint 1: -30 °C
- Setpoint 2: -50 °C

