



Setup details

Unistat® 910w & Radleys 10-litre reactor

Temperature range: -90...250 °C
 Cooling power: 5.2 kW @ 250...-20 °C
 4.7 kW @ -40 °C
 3.1 kW @ -60 °C
 Heating power: 6.0 kW
 Hoses: 2x1.5 m; M30x1.5 (#6386)
 HTF: DW-Therm (#6479)
 Reactor: 10-litre jacketed glass reactor
 Reactor content: 7.5 litre M90.055.03 (#6259)
 Stirrer speed: 200 rpm
 Control: process

Unistat® 910w

Cooling a Radleys 10-litre reactor to T_{min}

Requirement

The diagram illustrates the performance of a Unistat 910w undergoing two set-point changes, the second set-point is entered to find out the lowest temperature that the Radleys 10-litre reactor jacket and process can achieve in this set-up (" T_{min} ").

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 7.5 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The temperature profile is programmed, controlled and recorded with "SpyControl" software. The Unistat 910w is connected to a 10-litre glass reactor with a pair of M30x1.5 hoses.

For the first segment the process temperature reach -60 °C in approx. 100 minutes. Then the minimum process temperature achieved was -81 °C with a jacket temperature of -84 °C.

